

Small Animal



# Video- Bronchoscope

Professional bronchoscopes



# WIRELESS Video Bronchoscope

With integrated video camera  
and LED light

21 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com

**Dr. Fritz**<sup>®</sup>  
Wireless Endoscopes



# Video Bronchoscope

Professional bronchoscopes



- ✓ Integrated sender for wireless video transmission
- ✓ Easy handling, quick cleaning & disinfection
- ✓ Fully immersible
- ✓ Up/down controlled
- ✓ Working channel

| Item No.  | Length | ø     | Working channel    | Transmission |
|-----------|--------|-------|--------------------|--------------|
| TXBH2-36V | 36cm   | 2,9mm | No working channel | wireless     |
| TXBH4-54V | 54cm   | 3,9mm | 1,2mm              | wireless     |
| TXBH5-76V | 75cm   | 5,9mm | 1,8mm              | wireless     |

## Set TXBH4-54V

- Bronchoscope TXBH4-54
- Charger
- Cleaning brush
- Density tester
- Grasping forceps
- Storage case
- 7" Monitor receiver

## ► Several applications

see page 13



Wireless transmission





# Video Bronchoscope

## ■ Set TXBH4-54V

Bronchoscope  
OD: 3,9mm, WL: 54cm,  
Fully immersible  
(TXBH4-54)




Charger  
(TV27-NET)



Grasping forceps  
(F12-811)



Cleaning brush  
(F10-1107)



Density tester  
(EF-18)



7" Monitor receiver  
(TXM7-58R)



Case  
56x41x11cm  
(K-5041-11)

