


Small Animal



Compact Rhinoscope-/Vaginoscope

Classic endoscopic equipment for pets- and small animals



WIRELESS Compact Rhinoscope/ Vaginoscope

One of the smallest units with
working channel worldwide



Compact Rhinoscope/Vaginoscope

Classic endoscopic equipment for pets and small animals



- ✓ Battery replaceable
- ✓ Control display on the endoscope
- ✓ Integrated sender 5,8 GHz
- ✓ Integrated LED light
- ✓ Easy handling & quick cleaning
- ✓ Fully immersible



Control display



Battery replaceable

■ Set (TA20-3SET)

Cordless generation with exchangeable battery, in-built LED's, video camera, HD sender for video transmission & 1,5" control display.

- Compact Rhinoscope/Vaginoscope
- 2 spare batteries
- Charging station
- Charger
- Biopsy forceps
- 3 way stopcock
- Catheter, 1,0mm
- Cleaning brush
- Protection sheath
- NEW: enables circulation flushing
- 10" monitor receiver

▷ Several applications

see page 13

Wireless video transmission





Compact Rhinoscope/Vaginoscope

Classic endoscopic equipment for pets and small animals

■ Set (TA20-3SET)


Compact Rhinoscope-/Vaginoscope
OD: 3,4x2,6mm, 8,5 Charr
WL: 20cm
Working channel: 1,2mm
Luer connection, fully immersible



Charger
(97-FU14EU)

Charging station
(TVG-4)

Spare Li-Battery
(TVG-36)



Biopsy forceps
OD: 1,0 mm
WL: 35cm
(B10-351)



Flushing & protection cannula
4mm
enables 2 way suction & irrigation
(TV20-3SP)



Cleaning brush



3 way stopcock
(Y-LL-1273H)



Catheter
OD: 1,0 mm
WL: 100cm
(F13-102A)



10" Monitor receiver
incl. recorder & accu
(TXM10-58R)

