

Small Animal



Video Compact System 7

Endoscope for general ambulatory use, pets & exotic species

NEW
VCS7



WIRELESS Video Compact System 7

The next generation for your pet practice



Video Compact System 7

Endoscope for general ambulatory use, pets & exotic species

- ✓ Wireless video transmission
- ✓ Independent of video tower and light source
- ✓ Without cable
- ✓ Easy handling & quick cleaning
- ✓ Integrated LED light

NEW ⇒ Integrated working channel Ø 1,2mm
Ø 3,4mm & WL: 6cm

NEW ⇒ Chip-in-the-Tip-Technologie
robust and unbreakable



General use for small Animal consulting rooms

- Rhinoscopy
- Laryngoscopy
- Dental endoscopy
- Otoscopy
- Vaginoscopy
- Dermoscopy
- Best compliance with the owner

■ Set 98-VCS7-SE

- VCS7, Ø 3,4mm, WL: 6cm, working channel Ø1,2mm
- 10" monitor receiver
- Application sheaths for ears and oral examination in pets
- Optical grasping forceps
- Flexible grasping forceps Ø1mm, WL: 16cm
- Ear curette
- Charging station, wall mount
- Antifog, cleaning brushes
- Stainless steel stand for silicone sleeves including charging station
- Inflation pump for dental and gyn





Video Compact System 7

Endoscope for general ambulatory use, pets & exotic species

Das Set

VCS 7

Ø 3,4mm, WL:6cm

Working channel Ø 1,2mm

Completely soakable



Stainless steel stand for
silicone sleeves including
charging station
(981-VK-RA5)



Ear curette, oval, rounded end
OD:2mm WL:14cm
(B18-121)



Cats sheath
OD:4,5mm, WL:6,7cm
(98-VCS6-3)



Silicon sheath
OD:5mm, WL:7cm
(98-VCS6-6S)



Dental sheath
OD:8mm, WL:7,3cm
(98-VCS6-4)



Tongue spatula
OD:7,3mm, WL:8cm
(98-VCS6-5)



Laser sheath
OD:5,4x5,0mm, WL:6,8cm
working channel: 3,6mm
(98-VCS6-9)



Optical grasping forceps
(98-VCS6-8F)




Antifog fluid
(E-22-71)



Flexible grasping forceps
OD: 1mm, WL:16cm
(F10-160)



Cleaning brushes



(F15-181E)



(DE-5-17)

Inflation pump, small
4L/min., max. 0,2 bar,
incl. silicon tube
(H150-1-1)



10" Monitor receiver
incl. recorder & accu
(TXM10-58R)





Several applications

✓ Multi receiver capable

✓ Automatic radio connection, all types of endoscopes (broncho/arthroscopes/laparoscopes)



Receiver - monitor



TXL2-580M

5" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery
- With arm strap

TXM7-60R

7" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery

NEW: Time and labelling stamps for patient data

TXM10-58R

10" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery

TXM10-60R

10" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery

NEW: Time and labelling stamps for patient data



M-LCD-22PC

21" LCD PC

- With in-built receiver and recorder

TXL1-7-NO2

Receiver for Apple & Windows

- Incl. USB grabber



M-LCD-21SK

21" Video monitor HD

- With in-built receiver for scopes

TXL1-583

TV receiver for monitor

- Video-Out

TGR-40

Multi HD radio receiver

- HDMI out



Additional equipment


Video Compact System 7

Suction-/irrigation pump

For lens cleaning, inflation and suction.
For small and large animal practice.

H140-3S

- Suction & irrigation pump (H140-3)
- 1 li. suction bottle (H160-AB-43)
- 1 li. flushing bottle (H160-AB-7)
- Tubes & transformer
- Sterile filter (AR-10-922)




Optical suction & irrigation sheath

Suction & irrigation sheath
for VCS 6 & 7
(98-VCS6-11)



Flexible biopsy forceps

- Ø1mm, WL:15cm
- Premium quality
(B10-150)



Silicon tips

Silicon tips
Autoclavable
20 pcs.
(981-VK-S)



Suitable for silicon sheath
98-VCS6-6S

Working funnel

- Big (98-VCS6-7T)
- 9mm
- To take swab samples (gyn)
and large foreign body graspers

