


Bovine



Compact Video Vaginoscope

Professional video vaginoscope, wireless, for diagnostic use and treatment in dairy cattle



WIRELESS Compact Video Vaginoscope

**worldwide, unique and up to now
unmatched easy in practical use**



Control display



Battery replaceable

125 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com

Dr. Fritz[®]
Wireless ■ Endoscopes



Compact Video Vaginoscope

Professional video vaginoscope, wireless, for diagnostic use and treatment in dairy cattle

■ For practitioners and universities



- ✓ Battery replaceable
- ✓ Control display on the endoscope
- ✓ Integrated sender 5,8 GHz
- ✓ Integrated LED light
- ✓ Easy handling & quick cleaning
- ✓ Fully immersible

Set TQ52-28RD

Compact Video Vaginoscope & speculum for cattle. Cordless generation with exchangeable battery, built-in camera with 5,8 GHz radio transmission, LED light, 1,5" display. Both completely soakable when vaginoscope and battery a separated.

- Compact Video Vaginoscope
- Obturator, blunt
- 2 spare batteries
- Charging station
- Charger
- Reducer adapter
- 10" monitor receiver

■ Accessories to be used as transvaginal laparoscope

TQ52-65

Fornix perforator / cross trocar, 26x570mm. Compatible with TQ52-28RD.

- Trocar & blade guide (TQ52-65-1/2)
- Spare blades, 12 pcs./pack (TQ52-66KL)
- Spare part set, hexagonal wrench, spring, protection caps (TQ52-65ER)



TQ52-64S

Ova ligator, Ring ligator, 64cm. Compatible with TQ52-28RD.

- Main part with handle (TQ52-64G)
- Lifter (TQ52-64H)
- Ligatur rings (TQ52-R50)
- Applicator (TQ52-64Y)





Compact Video Vaginoscope

Professional video vaginoscope, wireless, for diagnostic use and treatment in dairy cattle

■ Complete Set for vaginoscopy and transvaginal laparoscopy and spaying



- ✓ Battery replaceable
- ✓ Control display on the endoscope
- ✓ Integrated sender 5,8 GHz
- ✓ Integrated LED light
- ✓ Easy handling & quick cleaning
- ✓ Fully immersible

Set TQ52-28RD1

TVC-Set, transvaginal spay/castration.
Cordless generation with exchangeable battery,
built-in camera with 5,8 GHz radio transmission, LED light,
1,5" display. Both completely soakable when vaginoscope
and battery a separated.

- Compact Video Vaginoscope
- Obturator, blunt
- 2 spare batteries
- Charging station
- Charger
- Reducer adapter
- Fornix perforator/ cross trocar
- 2 x ova ligator
- 50 ligatur rings
- 5" arm monitor

■ Spare parts

Scalpel blades

No.20 for trocar TQ52-65,
individuallypacked, 12pcs/pack
(TQ52-66KL)

Spare rings

Passend zu TQ52-64S,
production sterile,
FDA biocompatible, 100pcs/pack
(TQ52-R100)

Reducer adapter

For TQ52-28, 28/10mm, luer (fem),
silicon cap 8mm (234-211),
and sealing cap (234-600)
(TQ52-E)

Wireless HD transmisson





Individual parts TQ52 series

Compact Video Vaginoscope
OD: 28mm , WL: 52cm
fully immersible
(TQ52-28)



Charger
(97-FU14EU)



Charging station
(TVG-4)

Spare Li-Battery
(TVG-36)



Obturator, blunt
(TQ52-281)



Reducer adapter
28/10mm, luer (fem)
(TQ52-E)



Ova ligator, WL:64cm
(SET: TQ52-64S)



Fornix perforator/ cross trocar
(SET: TQ52-65)



Spare rings,
50pcs. (TQ52-R50)
100pcs. (TQ52-R100)



Scalpel blades, 12pcs.
(TQ52-66KL)



Case,
Case lid can be used as a
disinfection tray
Dimensions: 84x27x11cm
(Y-K-8427-11-1 /
Y-TQ52-SC)



5" Monitor receiver
incl. recorder & accu
(TXL2-580M)



10" Monitor receiver
incl. recorder & accu
(TXM10-58R)





Several applications

- ✓ Multi receiver capable

- ✓ Automatic radio connection, all types of endoscopes (broncho/arthroscopes/laparoscopes)



Receiver - monitor



TXM10-58R

10" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery

TXM7-60R

7" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery

TXL2-580M

5" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery
- With arm strap



M-LCD-22PC

21" LCD PC

- With in-built receiver and recorder

TXL1-7-NO2

Receiver for Apple & Windows

- Incl. USB grabber



M-LCD-21SK

21" Video monitor HD

- With in-built receiver for scopes

TXL1-583

TV receiver for monitor

- Video-Out

TGR-40

Multi HD radio receiver

- HDMI out