

Small Animal



Video- Bronchoscope

Professional bronchoscopes



WIRELESS Video Bronchoscope

With integrated video camera
and LED light

20 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com

Dr. Fritz[®]
Wireless Endoscopes



Video Bronchoscope

Professional bronchoscopes



- ✓ Integrated sender for wireless video transmission
- ✓ Easy handling, quick cleaning & disinfection
- ✓ Fully immersible
- ✓ Up/down controlled
- ✓ Working channel

Item No.	Length	ø	Working channel	Transmission
TXBH2-36V	36cm	2,9mm	No working channel	wireless
TXBH4-54V	54cm	3,9mm	1,2mm	wireless
TXBH5-76V	75cm	5,9mm	1,8mm	wireless

Set TXBH4-54V

- Bronchoscope TXBH4-54
- Charger
- Cleaning brush
- Density tester
- Grasping forceps
- Storage case
- 7" Monitor receiver

► Several applications

see page 13



Wireless transmission





Video Bronchoscope

■ Set TXBH4-54V

Bronchoscope
OD: 3,9mm, WL: 54cm,
Fully immersible
(TXBH4-54)




Charger
(TV27-NET)



Grasping forceps
(F12-811)



Cleaning brush
(F10-1107)



Density tester
(EF-18)



7" Monitor receiver
(TXM7-58R)



Case
56x41x11cm
(K-5041-11)

