



NEW!

Video hook grasper

The first video foreign body forceps for small and large patients

- For large item retrieval
- Rounded and self closing jaws
- Distal end curved
- Easy connectable to any video monitor
- Luer air connection




Video 3-hook grasper set, large

FK106-21S

- Ø 22mm, WL: 1,1m
- 3m video-/power cable & power adapter (110-230VAC;50/60Hz)
- Video cable with cinch
- Not soakable

Set consisting of:

- Video 3-hook grasper (FK106-21)
- Inflation pump (H150-1-1)
- 7" monitor with accu (TXM7-59A)




Video 2-hook grasper set, small

FK187-19S

- Ø 25x19mm, 65 Charr, WL: 87cm
- 3m video-/power cable & power adapter (110-230VAC;50/60Hz)
- Video cable with cinch
- Not soakable

Set consisting of:

- Video 2-hook grasper (FK187-19)
- Inflation pump (H150-1-1)
- 7" monitor with accu (TXM7-59A)



Video hook grasper also available individually

