



LED Dentalstick HD „Water Jet“

wireless, digital transmission



↖ Luer connection for manual and automatic flushing

The first dental endoscope with integrated irrigation

- Kabellos
- 5,8 GHz Digitalfunk
- HD Videoqualität
- 3 LEDS für die perfekte Ausleuchtung
- Integrierte Spülung
- Luer Anschluss
- 90° Seitblick und 140° Öffnungswinkel



↖ Opening for irrigation

■ Dentalstick professional HD complete set



TD-92HDS

- Dentalstick HD „Water Jet“ set
- Ø 12mm, NL:46cm, Wireless
 - Cordless high end video endoscope, in-built LED light and 5,8GHz sender rechargeable Li battery, resolution: PAL, HD 1280x720pxs
 - 140° field of vision
 - with integrated flushing
 - Luer connector (f)
- Set consisting of:**
- Dentscope HD
 - Charger
 - Short silicone tube with Luer connection on both sides (m/f)
 - 7" receivers monitor, DVR
 - Case

Schnurlose Videoübertragung

