


Equine



Endoscopic equipment

Nose, Larynx, Teeth



Endoscopy

for large Animals, Equine, Bovine & Camels.
Nose, Larynx, Teeth.





Dental - Large Animals

special instruments

LED illuminated mirror

for dental and pharyngeal examination in large animals

- With integrated LED light
- Robust
- Easy handling



Standard version D: 50mm

EB-95-5SET

LED Dental mirror Set large

D: 50 mm, WL:49 cm, complete with bag, anti-glare construction

EB-95-3SET

LED Dental mirror Set, small

D: 38 mm, WL:49 cm, complete with bag, anti-glare construction

Special Elevator

- Robust
- Easy cleaning



I-18-538

Special elevator long

D: 3mm, WL:20 cm, strong version




Flexibles LED Video Sinus Endoscope

Wireless video transmission for monitors or PC (<20m)

- Integrated video transmitter
- 2 way defecating
- 1,8 mm working channel
- Towerless application
- Fully immersible
- Easy and quick reprocessing

Video transmission



TXBH5-76V

Wireless video bronchoscope,

OD: 5,9mm, WL: 75cm

in-built 5,8 GHz sender for wireless video transmission and Li Ion battery resolution (400x400)

Set consisting of:

- Bronchoscope
- Charger
- Cleaning brush
- Density tester, case
- 7" receiver monitor





LED Dentalstick HD „Water Jet“

wireless, digital transmission



Luer connection for manual and automatic flushing

The first dental endoscope with integrated irrigation

- Wireless
- 5,8 GHz sender
- HD video quality
- 3 LEDS for perfect illumination
- Integrated flushing
- Luer connection
- 90° side view und 140° field of vision



Opening for irrigation


Dentalstick professional HD complete set



TD-92HDS

- Dentalstick HD „Water Jet“ set
- Ø 12mm, NL:46cm, Wireless
 - Cordless high end video endoscope, in-built LED light and 5,8GHz sender
 - rechargeable Li battery, resolution: PAL,HD 1280x720pxs
 - 140° field of vision
 - with integrated flushing
 - Luer connector (f)
- Set consisting of:**
- Dentscope HD
 - Charger
 - Short silicone tube with Luer connection on both sides (m/f)
 - 7" receivers monitor, DVR
 - Case
 - 7" Empfänger-Monitor mit Rekorder
 - Transportkoffer

Wireless transmission



Equine



LED-Dentalstick ECO

wireless, digital transmission

Finally back
in stock!



- Wireless
- 2 LEDs
- 90° side view und 120° field of vision
- With case

■ Dentalstick ECO complete set



TD-93ECS

Dentalstick ECO SET

- Ø 14mm, WL:46cm, wireless
- Resolution: PAL 768x756px
- Length: 70cm
- 120° field of vision

Set **consisting of:**

- Dentalstick ECO
- Charger
- 5" Monitor
- Case

Wireless transmission

