

Small Animal



Video- Gastroscope

Professional gastroscopes - connection for insufflation, irrigation and suction



WIRELESS Video Gastroscope

- NextGen4 -



23 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com

Dr. Fritz[®]
Wireless ■ Endoscopes



Video- Gastroscope

Professional gastroscopes - connection for insufflation, irrigation and suction

- ✓ Integrated sender for wireless video transmission
- ✓ Integrated video camera with LED light
- ✓ Integrated DVR
- ✓ Ergonomic body
- ✓ Easy handling, quick cleaning and disinfection


New Generation - NextGen4

- Integrated control display
- Integrated SD card slot
- Fully soakable
- Truly mobile

Item No.	Length	Ø	Working channel	Transmission
TG107-8	105cm	7,9mm	2,8mm	Wireless
TG136-9	135cm	9,8mm	2,8mm	Wireless
TG176-9	165cm	9,8mm	2,8mm	Wireless

Set TG136-9G

- Gastroscope TG136-9
- Charger
- Cleaning adaptors
- Cleaning brush
- Density tester
- Biopsy forceps
- 10" receiver-monitor, recorder
- Air/water pump
- Transportation case
- Optional: wall mount



Wireless HD transmission





■ Set TG136-9G

Gastroscope
OD: 9,8mm
WL: 1,35m
(TG136-9)




Charger
magnetic
(TV27-NET)



Cleaning adaptors
(TW132-24S)



5x Cleaning brush
(F28-280E)



Biopsy forceps
(B23-180E)



Density tester
(EF-23)



Air/water pump
(TW4-250SE)



10" Receiver-monitor,
recorder
(TXM10-58R)



Case for NextGen4
size: 56x41x11cm
(K-5641-10)





Several applications

- ✓ Multi receiver capable

- ✓ Automatic radio connection, all types of endoscopes (broncho/arthroscopes/laparoscopes)



Receiver - monitor



TXM10-58R

10" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery

TXM7-60R

7" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery

TXL2-580M

5" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery
- With arm strap



M-LCD-22PC

21" LCD PC

- With in-built receiver and recorder

TXL1-7-NO2

Receiver for Apple & Windows

- Incl. USB grabber



M-LCD-21SK

21" Video monitor HD

- With in-built receiver for scopes

TXL1-583

TV receiver for monitor

- Video-Out

TGR-40

Multi HD radio receiver

- HDMI out



Additional equipment


Video Gastroscope

Air/water pump

For lens cleaning and inflation

TW4-250SE

- Air/water pump (H150-1-1)
- Air/water twin tube (TW4-133)
- Fixation for bottle (TW4-261H)
- Bottle 250ml (TW4-250)




Suction-/irrigation pump

For lens cleaning, inflation and suction.
For small and large animal practice.

H140-3S

- Suction & irrigation pump (H140-3)
- 1 li. suction bottle (H160-AB-43)
- 1 li. flushing bottle (H160-AB-7)
- Tubes & transformer
- Sterile filter (AR-10-922)



Wall mount, inductive charging

For perfect storage

TGX-25S

- Wall mount
- Power supply





Additional equipment

Video Gastroscope

■ Catheter, injection needle, guttural pouch probe


F20-220

Teflon catheter
with luer-lock connector

- OD: 2 mm, WL: 235 cm
- autoclavable

F20-219

Spare catheter

- OD: 2 mm, WL: 10 m
- autoclavable

LLW-1220

Luer-lock connector

- Female OD: 2 mm
- Autoclavable


F22-16045


Guttural pouch probe

- OD: 2,2 mm, WL: 165 cm

F22-22046

Guttural pouch probe


- OD: 2,2 mm, WL: 230 cm


F231-220

Injection needle with teflon tube,
luer-look connector

- OD: 2,2 mm, WL: 235 cm,
Needle: 0,9x6 mm

■ Biopsy-grasping forceps, hooks, brush


B22-16011

Biopsy forceps with spike

- OD: 2,4 mm WL: 165 cm

B22-22011


Biopsy forceps with spike

- OD: 2,4 mm WL: 230 cm

B23-3601

Biopsy forceps with spike

- OD: 2,3 mm WL: 360 cm


F22-16041

Grasping forceps

- OD: 2,2 mm WL: 160 cm

F22-22041

Grasping forceps

- OD: 2,2 mm WL: 230 cm

F23-3604

Grasping forceps

- OD: 2,3 mm WL: 360 cm


B23-180E

Disposable biopsy forceps
steril

- OD: 2,3 mm WL: 180 cm

B23-230E

Disposable biopsy forceps
steril

- OD: 2,3 mm WL: 230 cm


F23-16042

Grasping forceps
„foreign body“

- OD: 2,3 mm WL: 160 cm

F23-22042

Grasping forceps
„foreign body“

- OD: 2,3 mm WL: 235 cm



Additional equipment

Video Gastroscope

FK23-1635

Basket size: 4,5 x 6cm

6 wires

- OD: 2,2mm WL: 160 cm



FK23-1637

Basket size: 5 x 10cm

6 wires

- OD: 2,3mm WL: 160 cm

FK23-2035

Basket size: 4,5 x 6cm

6 wires


- OD: 2,3mm WL: 230 cm

F22-16040

Grasping forceps

3 hooks

- OD: 2,2 mm WL: 165 cm



F22-22040

Grasping forceps

3 hooks

- OD: 2,2 mm WL: 230 cm

F28-1709

Cleaning brush

- OD: 2,6-3,5mm WL: 150 cm



F28-2208

Cleaning brush

- OD: 2,6-3,5mm WL: 230 cm



F28-3008

Cleaning brush

- OD: 2,6-3,5mm WL: 330 cm

F28-280E

Disposable cleaning brush

steril, 10pcs./pack

- OD: 1,8-2,8mm WL: 210 cm




For high frequency surgery

SL23-2310

Asymetric polypectomy snare complete with tube, without handle


- OD: 2,2mm, WL: 230cm wide: 25mm



SL23-1612

Symetric polypectomy snare complete with tube, without handle

- OD: 2,2mm, WL: 165cm wide: 35mm



SL23-2312

Symetric polypectomy snare complete with tube, without handle

- OD: 2,2mm, WL: 220cm wide: 35mm

HF22-1601

Suction-coagulation electrode monopolar, luer-lock connector

- OD: 2,2 mm, WL: 165 cm



HF-301-ER

HF- connecting cord

monopolar, com. to Erbe units 4mm

- Lenght: 3 m



HF22-2301

Suction-coagulation electrode monopolar, luer-lock connector

- OD: 2,2 mm, WL: 235 cm

HS11-1200

Handle for polypectom, HF snares

- standard connections





Additional equipment

Video Gastroscope


■ Brugman's foreign body catcher

for removing of foreign bodies
compatible to gastroscopes


FK-15-90 Brugman's foreign body catcher, short
90cm, OD:24mm, single use, 1 piece

FK-15-135 Brugman's foreign body catcher, long
135cm, OD:24mm, ID:15mm, 1 piece

compatible to gastroscopes
different manufacturers with OD <13mm



Gastroscope



* patent pending

- Attention: disposable product
return shipment not accepted.

X-46-00 Foreign body protection cap
for gastroscopes Ø 9-11 mm





Additional equipment

Video Gastroscope

■ For rigid and flexible endoscopes

DE-10-702 Dr. Fritz enzymatic cleaning flex,
1,5%, 1 litre, 7 enzymes, 0,5%,
for daily manual and automated cleaning



DE-10-300 Disinfection solution ,
for flexible and rigid endoscopes.
for manual and automated use,
1-3% dilutions, 5 litres




■ Option for cleaning of flexible endoscopes in the field

Disinfection-, transportation & storage basin

for flexible endoscopes

- For manual & semi-automatic cleaning
- Can also be used for cleaning machines
- Favourable
- For all endoscope types
- Can be integrated as a drawer
- With air vent
- Stackable



Stackable



DE-17-44S Basin for cleaning, transportation and
storage for flexible endoscopes
50x44x16 cm,

consisting of:
Basin (DE-17-44)
Cover, transparent (DE-17-41)



NEW!

Video hook grasper

The first video foreign body forceps for small and large patients

- For large item retrieval
- Rounded and self closing jaws
- Distal end curved
- Easy connectable to any video monitor
- Luer air connection




Video 3-hook grasper set, large

FK106-21S

- Ø 22mm, WL: 1,1m
- 3m video-/power cable & power adapter (110-230VAC;50/60Hz)
- Video cable with cinch
- Not soakable

Set consisting of:

- Video 3-hook grasper (FK106-21)
- Inflation pump (H150-1-1)
- 7" monitor with accu (TXM7-59A)




Video 2-hook grasper set, small

FK187-19S

- Ø 25x19mm, 65 Charr, WL: 87cm
- 3m video-/power cable & power adapter (110-230VAC;50/60Hz)
- Video cable with cinch
- Not soakable

Set consisting of:

- Video 2-hook grasper (FK187-19)
- Inflation pump (H150-1-1)
- 7" monitor with accu (TXM7-59A)



Video hook grasper also available individually

