

Small Animal



Compact Video-Vaginoscope TCI

Set universal for all dog sizes - transcervical insemination acc. to Blendingner



WIRELESS Compact Video-Vaginoscope

Special TCI system for transcervical
insemination acc. to Blendingner



Compact Video-Vaginoscope TCI

Set universal for all dog sizes - transcervical insemination acc. to Blendingner



- ✓ Battery replaceable
- ✓ Control display on the endoscope
- ✓ Integrated sender 5,8 GHz
- ✓ Integrated LED light
- ✓ Easy handling & quick cleaning
- ✓ Fully immersible

Working channel with
Luer connection

Control display

Battery replaceable

Compact Video-Vaginoscope & special TCI system for transcervical insemination, acc. to Blendingner, 36cm micro endoscope with integrated video camera and LED light and wireless video transmission.

■ Set (SET-GY-5VI)

- Compact Video-Vaginoscope
- 2 spare batteries
- Charging station
- Charger
- Sheath D:7mm
- Sheath D: 13mm
- Protection sheath
- 3 way stopcock
- Silicon cap
- Catheter, 1,6mm
- Cleaning brush
- 10" monitor receiver

▷ Several applications

see page 13

Wireless video transmission





Compact Video-Vaginoscope TCI

Set universal for all dog sizes - transcervical insemination acc. to Blendingner

Set

Compact Video-Vaginoskop (TA36-4-3)
OD: 4,9 mm x 3,6mm, 13 Charr
WL: 36cm
Working channel: 2,6mm
Battery replaceable



Charger
(97-FU14EU)



Charging station
(TVG-4)



Sheath for vaginoscopy
OD:7mm, 29cm
(TA36-45S)




Sheath for vaginoscopy
OD:13mm, 30cm
(TA36-49S)



3 way stopcock
(Y-LL-1273H)




Special catheter
OD: 1,6 mm, WL: 69cm
Angled tip, olive end
(F12-69)



10" Monitor receiver
incl. recorder & accu
(TXM10-58R)



Cleaning brush
OD: 2,0 mm
WL: 115cm
(F20-181E)






Accessories

Video Vaginoscope

Special catheter
OD: 1,3 mm, WL: 69cm
angled tip, olive end
(F11-69)

Special catheter
OD: 1,6 mm, WL: 69cm
angled tip, olive end
(F12-69)



Special catheter for insemination
OD: 1,6 mm, WL: 50cm
straight end, connection LL(f)
(F16-50)



Flexible biopsy forceps
OD: 2 mm, WL: 35cm
(B18-33)



Three way stopcock
plastic, disposable, non steril
Package of 10 pieces
(LL-1273H)



One way stopcock
plastic, disposable, non steril
Package of 10 pieces
(LL-1271H)



Sealing cap for working channel
opening 1mm,
Package of 5 pieces
(TR-60-1D)



Transportation case for vaginoscope
without picture
(K-3325-62)



Gynecology, Vaginoscopy & AI

Rigid endoscopes for medium and large dogs

■ Video Vaginoscope with sheath



TV27-5SET

Video Vaginoscope/Laparoscope
0°, OD: 5mm, WL: 27cm,
cordless generation with exchangeable battery,
with LED light, 5,8 Ghz

Consisting of:

- Compact Vaginoscope/Laparoscope (TA27-5)
- Charging station (TVG-4)
- 2x battery (TVG-36)
- Outer sheath 12mm (STV27-51)
working channel: 5,6mm, WL: 27cm
- Vaginoscope sheath (STV27-5V)
- Small air pump (H150-1-1)
- Power adapter (TV27-NET)
- 10" Receiver monitor (TXM10-58R)

Working channel



LED light

