

Equine



OP Equipment

EquiPorter

EquiPorter®

by F. Schlöderer & R. Fritz

www.equiporter.com



- the EquiPorter Pre OR safety sling -



110 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com

Dr. Fritz®
Wireless Endoscopes



OP Equipment

EquiPorter

■ Pre OR safety Sling**

A new development of a protection and handling sling system for equine general surgery

- To be mounted before surgery
- To take down patient during sedation and anesthesia
- To handle during surgery
- To assist during recovery/standing after surgery



EquiPorter®
by F. Schlegel & R. Fritz www.equiporter.com



4 point carrying frame



On-load operation



Belt elongation

- Easy mounting
- 4 point lifting/carrying system
- Individually adjustable
- Two sizes
- Washable

EP1-6170S

Pre OR safety belt, 170cm,
small, max.weight: 600kg
set incl.bag, washable (EP1-63)

EP1-6140S


Pre OR safety belt, 150cm,
small, max.weight: 500kg
set incl.bag, washable (EP1-63)

EP-R-63S

Carrying frame for OR safety belt
set with foam buffer,safety hooks, max. weight: 750kg
size: 101 x 40 x 15cm, 8,5kg



In combination with a 4 point carrying frame it will be
an ideal safety accessory for any
equine hospital with typical hoisting cranes



**Only for short time assisted lifting and weight support procedures (max.5-10min)

111 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com