# **Small Animal**



- Postrhinoscopy & Otoscopy
- Vaginoscopy
- Bronchoscopy & Gastroscopy
- Gynecology
- Arthroscopy
- Laparoscopy & Thoracoscopy



We are changing veterinary endoscopy



NEW VCS7 Video Compact System 7

The next generation for your pet practice





# **Video Compact System 7**

Endoscope for general ambulatory use, pets & exotic species

- √ Wireless video transmission
- √ Independent of video tower and light source
- √ Without cable
- √ Easy handling & quick cleaning
- √ Integrated LED light

Chip-in-the-Tip-Technologie robust and unbreakable



General use for small Animal consulting rooms

- Rhinoscopy
- Laryngoscopy
- Dental endoscopy
- Otoscopy
- Vaginoscopy
- Dermoscopy
- Best compliance with the owner

## Set 98-VCS7-SE

- VCS7, Ø 3,4mm, WL: 6cm, working channel Ø1,2mm
- 10" monitor receiver
- Application sheaths for ears and oral examination in pets
- Optical grasping forceps
- Flexible grasping forceps Ø1mm, WL: 16cm
- Ear curette
- Charging station, wall mount
- Antifog, cleaning brushes
- Stainless steel stand for silicone sleeves including charging station
- Inflation pump for dental and gyn







# **Video Compact System 7**

Endoscope for general ambulatory use, pets & exotic species

#### **Das Set**

VCS 7 Ø 3,4mm, WL:6cm Working channel Ø 1,2mm Completely soakable



Stainless steel stand for silicone sleeves including charging station (981-VK-RA5)



Ear curette, oval, rounded end OD:2mm WL:14cm (B18-121)

Cats sheath OD:4,5mm, WL:6,7cm (98-VCS6-3)



Silicon sheath OD:5mm, WL:7cm (98-VCS6-6S)



Dental sheath OD:8mm, WL:7,3cm (98-VCS6-4)



Tongue spatula OD:7,3mm, WL:8cm (98-VCS6-5)



Laser sheath OD:5,4x5,0mm, WL:6,8cm working channel: 3,6mm (98-VCS6-9)



Optical grasping forceps (98-VCS6-8F)



Antifog fluid (E-22-71)



Flexible grasping forceps OD: 1mm, WL:16cm (F10-160)



Cleaning brushes





Inflation pump, small 4L/min.,max.0,2 bar, incl. silicon tube (H150-1-1)



10" Monitor receiver incl. recorder & accu (TXM10-58R)



12 2025





# Several applications

√ Multi receiver capable



√ Automatic radio connection, all types of endoscopes (broncho/arthroscopes/laparoscopes)



**Receiver - monitor** 



#### TXL2-580M

5" Receiver-monitor, recorder

- mSD card slot, automatic scan, battery
- With arm strap

#### **TXM7-60R**

7" Receiver-monitor, recorder

mSD card slot, automatic scan, battery

NEW: Time and labelling stamps for patient data

#### TXM10-58R

10" Receiver-monitor, recorder

mSD card slot, automatic scan, battery

#### TXM10-60R

10" Receiver-monitor, recorder

mSD card slot, automatic scan, battery

NEW: Time and labelling stamps for patient data





#### M-LCD-22PC

21" LCD PC

· With in-built receiver and recorder

#### **TXL1-7-NO2**

Receiver for Apple & Windows

Incl. USB grabber





#### M-LCD-21SK

21" Video monitor HD

• With in-built receiver for scopes

#### TXL1-583

TV receiver for monitor

Video-Out

#### **TGR-40**

Multi HD radio receiver

HDMI out

13 2025





# Additional equipment

**Video Compact System 7** 


# **Suction-/irrigation pump**

For lens cleaning, inflation and suction. For small and large animal practice.

#### H140-3S

- Suction & irrigation pump (H140-3)
- 1 li. suction bottle (H160-AB-43)
- 1 li. flushing bottle (H160-AB-7)
- Tubes & transformer
- Sterile filter (AR-10-922)





## **Optical suction & irrigation sheath**

Suction & irrigation sheath for VCS 6 & 7 (98-VCS6-11)





## Flexible biopsy forceps

- Ø1mm, WL:15cm
- Premium quality (B10-150)





## Silicon tips

Silicon tips Autoclavable 20 pcs. (981-VK-S)





# **Working funnel**

- Big (98-VCS6-7T)
- 9mm
- To take swab samples (gyn) and large foreign body graspers







# Small Animal Middle ear diagnostic

Measurement of middle ear function in seconds for daily routine



- Examination and visual video representation of the external auditory canal with creation of a report in conjunction with the VCS
- Measurement of eardrum/middle ear function using a mini sound probe Confirmation of correct positions within seconds
- Comparative documentation, tympanogram, and photo output





#### The set for immediate use

TY25-SET-BA

- Basis unit with display
- Wall mount plate
- Power supply
- Foot switch
- Measurement hand piece
- Accessory kid with soft ear caps

⇒ Expandable with Video Compact System 6&7

For veterinary use only



# Compact Rhinoscope-/Vaginoscope

Classic endoscopic equipment for pets- and small animals



# Compact Rhinoscope/ Vaginoscope

One of the smallest units with working channel worldwide





# Compact Rhinoscope/Vaginoscope

Classic endoscopic equipment for pets and small animals



- √ Battery replaceable
- √ Control display on the endoscope
- √ Integrated sender 5,8 GHz
- ✓ Integrated LED light
- √ Easy handling & quick cleaning
- √ Fully immersible



Control display



Battery replaceable

# Set (TA20-3SET)

Cordless generation with exchangeable battery, in-built LED's, video camera, HD sender for video transmission & 1,5" control display.

- Compact Rhinoscope/Vaginoscope
- 2 spare batteries
- Charging station
- Charger
- Biopsy forceps
- 3 way stopcock
- Catheter, 1,0mm
- Cleaning brush
- Protection sheath

NEW: enables circulation flushing

• 10" monitor receiver



# > Several applications

see page 13




# Compact Rhinoscope/Vaginoscope

Classic endoscopic equipment for pets and small animals

### Set (TA20-3SET)





Biopsy forceps OD: 1,0 mm WL: 35cm (B10-351)



Flushing & protection cannula enables 2 way suction & irrigation (TV20-3SP)





3 way stopcock (Y-LL-1273H)



Catheter OD: 1,0 mm WL: 100cm (F13-102A)



10" Monitor receiver incl. recorder & accu (TXM10-58R)





# **Additional equipment**

**Compact Video Rhinoscope/Vaginoscope** 


## **Suction-/irrigation pump**

#### H140-3S

- Suction & irrigation pump (H140-3)
- 1 li. suction bottle (H160-AB-43)
- 1 li. flushing bottle (H160-AB-7)
- Tubes & transformer
- Sterile filter (AR-10-922)

For lens cleaning, inflation and suction. For small and large animal practice.





## **NEW! Post rhino video mirrow**

#### TA20-3P18

OD: 18mm, WL: 20cm With soft palate hook

Until now examination of retro pharynx was only possible for clinics and hospitals with appropriate endoscopic equipment. The newly developed video mirror system allows the general practiciones to perform retro pharyngoscopy. Grasping of foreign bodies or biopsies under visual control are possible.








# Post rhinoscopy LED mirror system

for cats and dogs



## Choanae mirror with LED light, independent of endoscopes



## Zubehör





**Antifog fluid** E-22-71

for endoscopes and mirrors, 50 ml

K-3325-61 Special transportation case

33x25x6 cm (without picture)

20 2025





**Professional bronchoscopes** 



With integrated video camera and LED light





# Video Bronchoscope

**Professional bronchoscopes** 



- ✓ Integrated sender for wireless video transmission
- ✓ Easy handling, quick cleaning & disinfection
- √ Fully immersible
- √ Up/down controled
- ✓ Working channel

Item No.	Length	Ø	Working channel	Transmission
TXBH2-36V	36cm	2,9mm	No working channel	wireless
TXBH4-54V	54cm	3,9mm	1,2mm	wireless
TXBH5-76V	75cm	5,9mm	1,8mm	wireless

## Set TXBH4-54V

- Bronchoscope TXBH4-54
- Charger
- Cleaning brush
- Density tester
- Grasping forceps
- Storage case
- 7" Monitor receiver

## > Several applications









#### Set TXBH4-54V







7" Monitor receiver (TXM7-58R)



Case 56x41x11cm (K-5041-11)







# Video- Gastroscope

Professional gastroscopes - connection for insufflation, irrigation and suction



# Video Gastroscope

- NextGen4 -

**Tips & tricks** Device care for flexible endoscopes









Professional gastroscopes - connection for insufflation, irrigation and suction



## **New Generation - NextGen4**

- Integrated control display
- Integrated SD card slot
- Fully soakable
- Truely mobile

### **Set TG136-9G**

- Gastroscope TG136-9
- Charger
- Cleaning adaptors
- Cleaning brush
- Density tester
- Biopsy forceps
- 10" receiver-monitor, recorder
- Air/water pump
- Transportation case Option

onal: wall mount	Air/water valve
Suction valve	/ / water valve
On/Off	Working channel
Record	U/D wheel
Brake U/D	
Blake 0/D	Brake L/R
SD card slot (DVR)	L/R wheel
Suction channel	\ Air/water coupling
	All / water coupling







VELV



#### **Set TG136-9G**



Charger magnetic (TV27-NET)





Cleaning adaptors (TW132-24S)



5x Cleaning brush (F28-280E)





Density tester (EF-23)



Air/water pump (TW4-250SE)



10" Receiver-monitor, recorder (TXM10-58R)



Case for NextGen4 size: 56x41x11cm (K-5641-10)







# Several applications

Multi receiver capable



Automatic radio connection, all types of endoscopes (broncho/arthroscopes/laparoscopes)



**Receiver - monitor** 



#### TXL2-580M

5" Receiver-monitor,recorder

- mSD card slot, automatic scan, battery
- With arm strap

#### **TXM7-60R**

7" Receiver-monitor, recorder

mSD card slot, automatic scan, battery

NEW: Time and labelling stamps for patient data

10" Receiver-monitor, recorder

mSD card slot, automatic scan, battery

#### TXM10-60R

10" Receiver-monitor, recorder

mSD card slot, automatic scan, battery

NEW: Time and labelling stamps for patient data



#### M-LCD-22PC

21" LCD PC

· With in-built receiver and recorder

#### **TXL1-7-NO2**

Receiver for Apple & Windows

Incl. USB grabber





#### M-LCD-21SK

21" Video monitor HD

With in-built receiver for scopes

#### TXL1-583

TV receiver for monitor

Video-Out

#### **TGR-40**

Multi HD radio receiver

HDMI out

27 2025





# Additional equipment

**Video Gastroscope** 


## Air/water pump

For lens cleaning and inflation

#### TW4-250SE

- Air/water pump (H150-1-1)
- Air/water twin tube (TW4-133)
- Fixation for bottle (TW4-261H)
- Bottle 250ml (TW4-250)






## **Suction-/irrigation pump**

For lens cleaning, inflation and suction. For small and large animal practice.

#### H140-3S

- Suction & irrigation pump (H140-3)
- 1 li. suction bottle (H160-AB-43)
- 1 li. flushing bottle (H160-AB-7)
- Tubes & transformer
- Sterile filter (AR-10-922)





#### Wall mount

Wall mount for inductive charging - always ready for use, perfectly stored and optimally drying.

#### **TGX-25S**

- Wall mount
- Power supply
- Mounting material







## Catheter, injection needle, guttural pouch probe



#### F20-220

Teflon catheter with luer-lock connector

- OD: 2 mm, WL: 235 cm
- autoclavable

#### F20-219

Spare catheter

- OD: 2 mm, WL: 10 m
- autoclavable

#### LLW-1220

Luer-lock connector

- Female OD: 2 mm
- Autoclavable



#### F22-16045

Guttural pouch probe

• OD: 2,2 mm, WL: 165 cm

#### F22-22046

Guttural pouch probe

• OD: 2,2 mm, WL: 230 cm



#### F231-220

Injection needle with teflon tube, luer- look connector


• OD: 2,2 mm, WL: 235 cm, Needle: 0,9x6 mm

## Biopsy-grasping forceps, hooks, brush

#### B22-16011

Biopsy forceps with spike

OD: 2,4 mm WL: 165 cm



#### B22-22011

Biopsy forceps with spike

OD: 2,4 mm WL: 230 cm

#### B23-3601


Biopsy forceps with spike

• OD: 2,3 mm WL: 360 cm

#### F22-16041

Grasping forceps

• OD: 2,2 mm WL: 160 cm



#### F22-22041

Grasping forceps

OD: 2,2 mm WL: 230 cm

#### F23-3604

Grasping forceps

• OD: 2,3 mm WL: 360 cm



#### B23-180E

Disposable biopsy forceps

• OD: 2,3 mm WL: 180 cm



Disposable biopsy forceps steril

• OD: 2,3 mm WL: 230 cm



#### F23-16042

Grasping forceps "foreign body"

OD: 2,3 mm WL: 160 cm



Grasping forceps "foreign body"

• OD: 2,3 mm WL: 235 cm





# Additional equipment

#### **Video Gastroscope**

#### FK23-1635

Basket size: 4,5 x 6cm 6 wires

• OD: 2,2mm WL: 160 cm



#### FK23-1637

Basket size: 5 x 10cm 6 wires

OD: 2,3mm WL: 160 cm

#### FK23-2035

Basket size: 4,5 x 6cm

6 wires

OD: 2,3mm WL: 230 cm



#### F22-16040

Grasping forceps 3 hooks

• OD: 2,2 mm WL: 165 cm



#### F22-22040

Grasping forceps 3 hooks

• OD: 2,2 mm WL: 230 cm



#### F28-1709

Cleaning brush


• OD: 2,6-3,5mm WL: 150 cm



# F28-280E

Disposable cleaning brush steril, 10pcs./pack

• OD: 1,8-2,8mm WL: 210 cm



#### F28-3008

F28-2208


Cleaning brush

Cleaning brush

• OD: 2,6-3,5mm WL: 330 cm



# For high frequency surgery



#### SL23-2310

Asymetric polypectomy snare complete with tube, without handle


OD: 2,2mm, WL: 230cm wide: 25mm



#### SL23-1612

Symetric polypectomy snare complete with tube, without handle

• OD: 2,2mm, WL: 165cm wide: 35mm



#### SL23-2312

Symetric polypectomy snare complete with tube, without handle

• OD: 2,2mm, WL: 220cm wide: 35mm



# HF22-1601

Suction-coagulation electrode monopolar, luer-lock connctor

• OD: 2,2 mm, WL: 165 cm



#### HF-301-ER

HF- connecting cord monopolar, com. to Erbe units 4mm

Lenght: 3 m



#### HF22-2301

Suction-coagulation electrode monopolar, luer-lock connctor

• OD: 2,2 mm, WL: 235 cm



#### HS11-1200

Handle for polypectom, HF snares

standard connections



30 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com



**Video Gastroscope** 


## Brugman's foreign body catcher

for removing of foreign bodies compatible to gastroscopes

FK-15-90 Brugman's foreign body catcher, short 90cm, OD:24mm, single use, 1 piece

FK-15-135 Brugman's foreign body catcher, long 135cm, OD:24cm, ID:15mm, 1 piece

> compatible to gastroscopes different manufacturers with OD <13mm







- \* patent pending
- Attention: disposable product return shipment not accepted.

X-46-00 Foreign body protection cap for gastroscopes Ø 9-11 mm







# Additional equipment

**Video Gastroscope** 

## For rigid and flexible endoscopes

#### DE-10-702 Dr. Fritz enzymatic cleaning flex,

1,5%, 1 litre, 7 enzymes, 0,5%, for daily manual and automated cleaning



#### DE-10-300 Disinfection solution,

for flexible and rigid endoscopes. for manual and automated use. 1-3% dilutions, 5 litres



## Option for cleaning of flexible endoscopes in the field

## Disinfection-, transortation & storage basin

for flexible endoscopes



- Can also be used for cleaning machines
- Favourable
- For all endoscope types
- Can be integrated as a drawer
- With air vent
- Stackable





Stackable

**DE-17-44S** 

Basin for cleaning, transportation and storage for flexible endoscopes 50x44x16 cm,

consisting of: Basin (DE-17-44)

Cover, transparent (DE-17-41)







# Video hook grasper

The first video foreign body forceps for small and large patients

- For large item retrieval
- Rounded and self closing jaws
- Distal end curved

NEW!

- Easy connectable to any video monitor
- Luer air connection







#### FK106-21S

- Ø 22mm, WL: 1,1m
- 3m video-/power cabel & power adapter (110-230VAC;50/60Hz)
  - -Video cable with cinch
- Not soakable

#### Set consisting of:

- Video 3-hook grasper (FK106-21)
- Inflation pump (H150-1-1)
- 7" monitor with accu (TXM7-59A)



#### FK187-19S

- Ø 25x19mm, 65 Charr, WL: 87cm
- 3m video-/power cabel & power adapter (110-230VAC;50/60Hz)
  - -Video cable with cinch
- Not soakable


#### Set consisting of:

- Video 2-hook grasper (FK187-19)
- Inflation pump (H150-1-1)
- 7" monitor with accu (TXM7-59A)





### Video hook grasper also available individually







33 2025



# **Compact Video-Vaginoscope TCI**

Set universal for all dog sizes - transcervical insemination acc. to Blendinger



# Compact Video-Vaginoscope

**Special TCI system for transcervical** insemination acc. to Blendinger



# **Compact Video-Vaginoscope TCI**

Small Animal

Set universal for all dog sizes - transcervical insemination acc. to Blendinger



Compact Video-Vaginoscope & special TCI system for transcervical insemination, acc. to Blendinger, 36cm micro endoscope with integrated video camera and LED light and wireless video transmission.

# Set (SET-GY-5VI)

- Compact Video-Vaginoscope
- 2 spare batteries
- Charging station
- Charger
- Sheath D:7mm
- Sheath D: 13mm
- Protection sheath
- 3 way stopcock
- Silicon cap
- · Catheter, 1,6mm
- Cleaning brush
- 10" monitor receiver

## > Several applications

see page 13










# Compact Video-Vaginoscope TCI

Set universal for all dog sizes - transcervical insemination acc. to Blendinger


#### Set





Sheath for vaginoscopy OD:7mm, 29cm (TA36-45S)

Sheath for vaginoscopy OD:13mm, 30cm (TA36-49S)





Special catheter OD: 1,6 mm, WL: 69cm Angled tip, olive end (F12-69)



10" Monitor receiver incl. recorder & accu (TXM10-58R)



Cleaning brush OD: 2,0 mm WL: 115cm (F20-181E)




Video Vaginoscope

Special catheter OD: 1,3 mm, WL: 69cm angled tip, olive end (F11-69)

Special catheter OD: 1,6 mm, WL: 69cm angled tip, olive end

(F12-69)



Special catheter for insemination OD: 1,6 mm, WL: 50cm straight end, connection LL(f) (F16-50)



Flexible biopsy forceps OD: 2 mm, WL: 35cm (B18-33)



Three way stopcock plastic, disposable, non steril Package of 10 pieces (LL-1273H)



One way stopcock plastic, disposable, non steril Package of 10 pieces (LL-1271H)



Sealing cap for working channel opening 1mm, Package of 5 pieces (TR-60-1D)



Transportation case for vaginoscope without picture (K-3325-62)





# Gynecology, Vaginoscopy & Al

Rigid endoscopes for medium and large dogs

## **Video Vaginoscope with sheath**




**TV27-5SET** 

Video Vaginoscope/Laparoscope 0°, OD: 5mm, WL: 27cm, cordless generation with exchangeable battery, with LED light, 5,8 Ghz

#### Consisting of:

- Compact Vaginoscope/Laparoscope (TA27-5)
- Charging station (TVG-4)
- 2x battery (TVG-36)
- Outer sheath 12mm (STV27-51) working channel: 5,6mm, WL: 27cm
- Vaginoscope sheath (STV27-5V)
- Small air pump (H150-1-1)
- Power adapter (TV27-NET)
- 10" Receiver monitor (TXM10-58R)





LED light





**Endoscopes and instruments for minimally invasive surgery** 



# Video-Laparoscope

for small animals





Telescopes for diagnostic and surgery

## Video Laparoscope, wireless




**TV27-5SET** 

Video Laparoscope/ Vaginoscope 0°, OD: 5mm, WL: 27cm, cordless generation with exchangeable battery, with LED light, 5,8 Ghz

#### Consisting of:

- Compact Vaginoscope/Laparoscope (TA27-5)
- Charging station (TVG-4)
- 2x battery (TVG-36)
- Outer sheath 12mm (STV27-51) working channel: 5,6mm, WL: 27cm
- Vaginoscope sheath (STV27-5V)
- Small air pump (H150-1-1)
- Power adapter (TV27-NET)
- 10" Receiver monitor (TXM10-58R)





LED light







Trocar cannulas with magnetic valve, safety trocars

Trocar cannulas are designed with different types of valves to reduce or to avoid CO2 loss when changing instruments and telescope during laparoscopy and thoracoscopy. Our reusable and patented magnetic valve trocar system has several advantages compared to other valve types. First of all they are light in weight and because of the valve construction, the insertion of instruments and scopes will pass smoothly through with out scratching or touching the front lens like in a ball system or in a so called automatic valve system. In addition there is nearly no gas loss even when changing instruments. There are no mechanical parts to clean inside except the magnetic flap itself which is a great help to your staff.

## Inner diameter 5,5mm

for standard 5mm instruments

231-213V Cannula only with magnetic valve, ID: 5,5 mm,

7 cm, stainless steel, stopcock, 40 g

Cannula only with magnetic valve, ID: 5,5 mm, 231-212V

7 cm, 30 g (without picture)

Safety trocar only "Cavi Safe" 231-137V

cork screw handles, compatible with cannulas

ID: 5,5 mm, WL: 7 cm, fully dismantable

231-105V Trocar only with sharp tip, compatible with

Canully only ID: 5,5 mm, WL: 7 cm (without picture)







Tissue pressure exposes the sharp tip

## Inner diameter 10mm

for telescopes and 10,, instruments

Cannula only with magnetic valve, 230-212V

ID:10 mm,WL: 7 cm, inflation stopcock, 55 g

Cannula only with magnetic valve, 230-213V

ID: 10 mm, WL: 7 cm, 45 g, (without picture)

230-137V Saftey trocar only with cork screw handles,

compatible to cannulas ID: 10 mm, WL: 7 cm,

fully dismantable

230-105V Trocar only with pyramidal tip for cannulas

ID: 10 mm, NL: 7 cm

## **Inner diameter 13mm**

for telescopes and 13mm instruments

637-232 Cannula only with magnetic valve,

ID:13 mm, WL: 10 cm, inflation stopcock, Trinox

Canulla only with pyramidal tip for cannulas 237-173

ID: 13 mm, WL: 10 cm





237-173





**Accessoriers, Trocar disposable** 



Canula only Magnetic valve flap Sealing caps		
231-211 Sealing caps for Ø 5,5mm, blue Package of 10 pieces	631-901  Magnetic valve body, Ø 5mm  With Magnetic Head-/valve flap & sealing caps	
234-211 Sealing caps for Ø 8mm, white Package of 10 pieces	634-901  Magnetic valve body, Ø 8mm  With Magnetic Head-/valve flap & sealing caps	
230-211 Sealing caps for Ø 10mm, red Package of 10 pieces	630-901  Magnetic valve body, Ø 10mm  With Magnetic Head-/valve flap & sealing caps	
231-215 Magnetic valve flap for Ø 5,5mm cannulas Package of 2 pieces	630-902 Spare stopcock Complete with spring nut	
230-215 Magnetic valve flap for Ø 8/10/11/12/13 mm cannulas Package of 2 pieces	630-903 Spare spring nut	
236-217 Sealing ring for Trinox trocar cannulas for Ø 3,5mm & 5,5mm Package of 5 pieces	230-217 Sealing ring for Trinox trocar cannulas for Ø 8, 10, 11, 12, 13mm Package of 5 pieces	
230-505 Snap reducer for cannulas from Ø 10mm to 5,5mm	237-510 Snap reducer for cannulas from Ø 13mm to 10mm	

# Trocar disposable 5,5mm & 12,5mm



TR-59-54

Trocar disposable, ID: 5,5mm, WL: 7 cm, with spike, LL(f), inflation stopcock,

package of 8 pieces

TR-60-121 Trocar disposable, ID: 12,5mm, WL: 10cm,

with optical spike, no blade, 12-8/-5 mm, LL(f) inflation stopcock,

package of 6 pieces

42 2025





**Palpation and Puncture** 

## **Laparoscopy first puncture needle**

L-19-100 Veress needle

OD: 2,1 mm, 10 cm, LL(f)





Injection cannula

Injection system and probe

L-210 Injection cannula

thin, 1x26 mm, 2 pieces

L-218 Injection cannula

large, 1,9x26 mm, 2 pieces

L-218R Injection cannula

rounded, 1,9x28 mm, 2 pieces

L-33-210 Injection system and probe

5 mm, 33 cm, protected interchangable needle,

consisting of:

Inside tube, LL, connector L-33

L-33-2 Outer tube/probe

L-210 Needle, 1 piece (1x26mm)

# Mini instrument for two portal technique in laparosopy, 2,4mm

L-24-23 Needle with mini grasping forceps

**O**D: 2,4 mm, WL: 23 cm





43 2025



Standard instruments, 5mm

We are currently optimizing our 5 mm instruments for you. Individual components may be available upon request — feel free to contact us!

#### **Plastic handles**

- With or without lock
- HF connector
- Colour coding for scissors


#### H 01/H 05



- Long life span
- **■** Compatible to 5mm trocar cannulas

#### **Tubular sheaths**

- 5 mm Ø
- Different working lengths
- Rotatable through 360° oder feststehend
- Insulated



#### Complete instruments, fully dismantable

#### **Scissors**

531.04H05V

Hook scissors, 5mm, 33cm, single action jaws, insulated, rotatable, dismantable, high quality, consisting of: jaws insert (53M04), sheath (53 R 1) handle (H05), color code:white



## **Biopsy forceps**

531.61H01V

Biopsy forceps, 5mm, 33cm, double action jaws, dismantable, consisting of: jaws insert (53M60), sheath (53R1) handle (H01), color code:black







## Laparoscopy / Thoracoscopy

Standard instruments, 5mm

We are currently optimizing our 5 mm instruments for you. Individual components may be available upon request — feel free to contact us!

#### **Universal dissection forceps**

#### 531.40H01V

Dissection forceps (Maryland) Ø 5 mm, WL: 33 cm

curved, insulated, rotatable, dismantable, consitsting of:

jaws insert (53M40), sheath (53R1) handle (H01) colour code: black



#### Dissecting forceps, strong Ø 5mm, WL: 33cm

curved, for use with 8 & 10mm trocars consitsting of: jaws insert (53M37), sheath (53R1) handle (H01) colour code: black



#### 531.22H01V

Grasping forceps, universal Ø 5 mm, WL: 33 cm

insulated, rotatable, zerlegbar, dismantable

consitsting of:

jaws insert (53M22), sheath (53R1) handle (H01) colour code: black



#### Bowel grasping forceps Ø 5 mm, WL: 33 cm

insulated, rotatable, zerlegbar, dismantable

consitsting of:

jaws insert (53M36), sheath (53R1) handle (H01) colour code: black



#### Grasping forceps insert Ø 5 mm, WL: 33 cm

front hooks consitsting of:

jaws insert (53M24), sheath (53R1) handle (H01) colour code: black



#### Atraumatic grasper Taccun, Ø 5 mm, WL: 33 cm

insulated, rotatable, zerlegbar, dismantable

consitsting of:

jaws insert (53M25), sheath (53R1) handle (H01) colour code: black



#### DeBakey grasping forceps Ø 5mm, WL: 33cm

consitsting of:

jaws insert (53M35), sheath (53R1) handle (H01) colour code: black





## Laparoscopy / Thoracoscopy

**Palpation** 

#### **Knot tier and cutter**



Easy handling, easy to tighten, additional cutting


L-19-534

Knot tier, standard version, OD: 5 mm, WL: 29 cm, easy filling, suture size up to 1metric

#### **Needle holders**

Axial handle, ergonomic, stable, O D: 5mm

- Easy Clean System
- Needle self positioning







### **Special devices - High Frequency Surgery**

High frequency unit, High speed, Typ Atom

High frequency unit, High speed - with vessel sealing technique. Parallel usage of all outputs.

- Adapter for instruments other manufacturers on request
- uncomplicated use
- up to date safety regulation
- vessel sealing with confirmation signal
- compatible to all instruments on the market
- program via Touch Screen
- SDS self detecting system



HF-401-TH HF/Radio frequency unit, High End Thermo Stapling (vessel sealing), touch-screen programmable, latest mono & bipolar technique SDS - self detecting system incl. high speed vessel sealing delivery without accessories, 400W, 230VAC, 50/60Hz

_				
Sn	ecif	いつさ	ION	יסו
$\mathbf{v}$	CUII	ıvaı	ıvı	ıo.

opeomedationer.			
Monopolar	Cutting Coagulation Soft/Forced	400 W 200 W + / +	
Bipolar	Cutting Coagulation Auto Start/Stop	400 W 200 W + / +	
Bipolar	Thermo-Stapling Thermo sealing Auto Stop	+ 300 W + +	
Weight	7 kg		
max. power		420 W	



### **Special devices - Inflator**

CO<sub>3</sub>Inflator, electronically monitored

- small and large animal
- reusable tubing
- bacterial filter
- economical



741-025V

CO<sub>2</sub> Inflator, electronically monitored, max. 20 litre, complete without bottle automatic pressure reduction system, DIN gas connection type, 115-230 VAC, 50-60Hz, DIN IEC 601/1, DIN VDE 871

#### Accessories





741-009-1 Spare Silicon tube, 3m, for CO, inflation,

incl. connector, autoclavable

741-009-FF Sterile filter system for CO<sub>2</sub> -inflator

disposable, box with 50 pcs. medical connection D:22mm

741-002 CO, high pressure tube, 1m

DIN connector, screw both sides

741-009V-PN CO, high pressure tube, 0,8m

compatible to PIN Index bottle (without picture)

Sterile product return shipment not accepted





### Laparoscopy/Thoracoscopy

Complete set with video laparoscope - wireless

#### Basic set without camera system and light source

#### SET-K-LA-B1

Laparoscopy set, small animals, conventional for stndard 2 portal technique complete, but without video camera and light source

#### consisting of:

- Wide angle telescope 30°/5mm/28cm (T281-5030)
- Safety trocar complete, 5,5mm (231-213V & 231-137V)
- Safety trocar complete, 10mm (230-212V & 230-137V)
- Trocar reducer 10/5 (230-505)
- Veress needle, 10cm (L-19-100)
- Grasping forceps, 5mm
- Metzenbaum scissors, 5mm
- Preparation forceps, 5mm
- CO2 Inflator, electronically monitored 20l (741-025V)



#### **Option:**

#### SET-K-LA-B2

same as SET-K-LA/VA1 + HF vessel sealing, type Atom incl. sealing MIC instrument 5mm

#### Complete set with video Laparoscope - wireless

#### SET-K-LA/VA1

Universal Laparo-/Vaginoscope SET, for all laparoscopic aproaches in small animals, for 1 & 2/3 portal technique and transvaginal laparoscopy

#### consisting of:

- LED Video Laparo. 0°/5mm/28cm (TA27-5)
- Working sheath 12mm (STV27-51)
- Safety trocar complete, 5,5mm (231-213V & 231-137V)
- Trocar disposable (TR-60-121)
- Veress needle, 10cm (L-19-100)
- Needle with mini grasping forceps (L-24-23)
- Grasping forceps, 5mm
- Metzenbaum scissors, 5mm
- Preparation forceps, 5mm
- CO2 Inflator, electronically monitored 20I (741-025V)
- 22" receiver/PC/recorder



#### **Option:** SET-K-LA/VA2

same as SET-K-LA/VA1 + HF vessel sealing, type Atom incl. sealing MIC instrument 5mm






Puncture System acc. to Dr.Lehmann for the canine shoulder joint

#### Innovation in small animal arthroscopy

Developed an tested at the Veterinary Faculty of Vienna, Austria

During recent years, as arthroscopy in small animals has developed, researchers realized the necessity of a specific puncture system to create an instrument portal in the larger joints. Due to the dense muscle mass in shoulder joints, creating a puncture for a working cannula is often quite difficult. If a clean puncture allowing access directly to the joint is not achieved the first time, the risk of a subsynovial oedema increases greatly, thus decreasing the time frame to complete the surgery with each failed puncture. The shoulder puncture system has been specifically designed to help eliminate this risk. The system allows you to align your arthroscope accurately in the joint space. Simply attach the puncture system to the arthroscopic sheaths and maintain your arthroscopic vision in a parallel direction to the shoulder joint rim. At this point you will find the working triangulation portal in an ideal position. After a puncture is created with the desired trocar, simply replace the trocar with a working cannula and start working with your instruments.





#### Irrigation and Working Cannulas, Trocar Changing System acc. to VAN BREE

The trocar changing system acc. to VAN BREE allows the specific exchange of different working cannulas. Beginning with the smallest diameter, the surgeon can always change to a bigger diameter by using the changing rods. According to the different joints a so called "elbow set" (working lengh 3,5-5 cm) and a "shoulder set"(7-7,5 cm) is differentiated. In addition it is possible to change the working portal with the telescope portal using the changing rod. The color code system allows easy size differentiation.

#### Set I-17-755

Working cannula set "shoulder and elbow" acc. to van BREE consisting of:

- Set I-17-500
- Set I-17-700





#### Set I-17-500

Working cannula set "elbow" acc. to van Bree, consisting of:

- I-17-524

Working cannula incl. sharp obturator ID: 2,3mm WL 5,5cm, gold

Working cannula ID: 2,9mm WL 4,8cm, metal

- I-17-535

Working cannula ID: 3,5mm WL 3,5cm, schwarz

Changing rod 2,3mm round / pointed, gold



I-17-024

#### Set I-17-700

Working cannula set "shoulder" acc. to van Bree, consisting of:

- I-17-728

Working cannula incl. sharp obturator ID: 2,9mm WL 7,5cm, metal

Working cannula ID: 3,5mm WL 6,8cm, black

- I-17-028

Changing rod 2,9mm round / pointed, metal



I-17-028

I-17-535-G

Working cannula with thread ID:3,5mm, WL:3cm (can replace I-17-535)



I-17-002

Sealing caps for small trocar cannulas opening: 2mm 5pcs.

I-17-735-G

Working cannula with thread ID:3,5mm, WL:6,8cm (can replace I-17-735)



51 2025





Hand instruments with colour coding are easy to recognize and perfectly to use in small and in large joints according to the patient sizes.

- **Easy handling**
- **Ergonomic design**
- For universal use
- Suitable for Dr. Frizu working cannulas



- AR-17-03S Bajonett knife 2,7mm, silicon handle, autoclavable



- AR-17-05S Banana knife, 2,7mm, silicon handle, autoclavable



- AR-17-06S Currette, small 2,7mm, silicon handle, autoclavable



- AR-17-07S Ring currette, 2,7mm, silicon handle, autoclavable



- AR-17-09S Hook probe, 2,2mm, silicon handle, autoclavable



- AR-17-11S Smilie knife, 2,7mm, silicon handle, autoclavable



- AR-17-12S Mikro picking knife, acc. to Dr. Beale, 2,7mm, silicon handle, autoclavable







#### **Hand instruments**

#### Elevator

AR-17-13

Special elevator, 2,7mm, with handle, metal



#### Milling Drill acc. to VAN RYSSEN & VAN BREE

The latest development allows cartilage and bone curretage even through the smallest working cannula

I-17-822

Milling drill for arthroscopy, acc.to van RYSSEN D: 2,2mm, WL: 12cm with rotating handpiece, long lasting drill steel, dismantable Constisting of: handpiece (I-17-820) and milling drill (I-17-821) Changeable drill

Rotating handpiece

### **Instruments for high frequency**

Mono-/Unipolar HF/RF cutting electrodes, the electric knife

HF-200-055

High frequency connecting cord comp. to arthroscopic cutting electrodes, handles length: 5 m, for units with 4 mm connector



HF-200-45

High frequency cutting electrode for arthroscopy, 45°, fine tip 4 mm connector, 12 cm



HF-200-90

High frequency cutting electrode for arthroscopy, 90°, heavy version 4 mm connector, 12 cm





**Punches, Biopsy Forceps, Rongeurs and Grasping Forceps** 

#### **2mm Instruments**

- AS-17-521

Mini scissors - high performance -D: 2mm, WL:11cm



#### 2,2/2,7 mm Instruments

AT-19-231

Mini-Rongeur - high performance -D: 2,2/2,7mm, WL:11cm Overstraining protection, cleaning connector



#### 2,7 mm Instruments

- AS-17-632

Universal rongeur and grasping forceps D: 2,7 mm, WL:12cm

- with overstraining protection -



#### - AS-17-641

Rongeur and biopsy forceps, D: 2,7 mm, WL: 12cm

- with overstraining protection -



- AS-17-672

Arthroscopic punch, D: 2,7 mm, WL:12cm

- with overstraining protection -



- AS-17-621

Arthroscopic punchh, D: 2,7 mm, NL:12cm



- AS-17-662 - mit Sperre -

Alligator grasping forceps,

D: 2,7 mm, WL:12cm

small grip



- Overstraining protected
- Wasp waisted design
- Long life span



54 2025





**Punches, Biopsy Forceps, Rongeurs and Grasping Forceps** 

#### 3,4 mm Instruments

- AT-19-343

Maxi-Rongeur -high performance-D:3,4 mm,WL:12 cm, with overstraining protection, cleaning connector



- AS-17-762

Rongeur, universal & grasping D:3,4mm,WL:12cm, with overstraining protection




#### 3,5 mm Instruments

- AY-19-463

Arthro grasper, grooves, D:3,5 mm, WL:10 cm, conical pointed, with maximum power U-shaped construction







- Fine and extremely strong
- **Ergonomically adapted with lock**

- AY-19-451 Arthro punch, straight, D:3,5 mm, W L:10 cm conical pointed, with maximum power cut

U-shaped construction





### 4,0 mm Instruments

AY-19-452

Arthro punch, bent upwards, D:4 mm, WL:10cm conical pointed, with maximum power cut U-shaped construction










#### Manual fluid supply system

One of the cheapest options for liquid supply



#### **Pressure Infusion Cuffs**



> For high flow fluid supply, please use big infusion tubes

#### AR-08-100N Pressure infusion cuff

for 500-1000ml with manometer and hand pump reusable high quality manufacturing

#### Roller pump for fluid supply

Large and small animal arthroscopy



- Quick and easy fluid flow
- Simple speed control
- No disposable accessories

#### ARR-20

#### Roller pump

Noise level: max. 50 dBa, very low Pump type: 3-roller peristaltic Flow parameters: up to 525ml/min.

Pneumatic foot switch

Weight: 4,4kg

Size: 190 x 230 x 330mm

#### **Accessories**



I-18-518 **Egress cannula** 

> with sharp trocar stopcock, ID: 2,7mm, WL: 7cm





#### Connection: Fluid bag to arthroscope

#### AR-10-948 Silicon tube

D: 6x2 mm, WL:3 m, autoclavable, Luer-Lock (m), needle f. sterile bags/bottles,

large inner diameter

X-7-20-10 Silicon spare tube

Silicon spare tube, 100cm, large

inner diameter, size: 6x2mm, color: natur

AR-10-9 Needle for infusion bottles/bags

ID: 4 mm, compatible to silicon tubes 3,5-5 mm

LLM-1258 Luer-Lock connector, male

10 mm


56 2025







#### Cordless single use shaver system for arthroscopy



#### **Starter set (SFR-205SET)**

Set consistinf of:

- Central unit motor (SFR-206M)
- Charging station (SFR-205B)
- Burr shaver Ø 4,5mm (SFR-104C)



#### SFR-206M Central unit - motor

The reusable central unit, the heart of the system, is equipped with state-of-the-art electronics and a rechargeable battery. It provides the surgeon with both standard and advanced features: a powerful brushless motor ensures optimal torque for all arthroscopic procedures, while the durable battery supports even the longest surgeries.

- 6000 rpm
- Approx. 40min. runnig time



#### SFR-205B Charging station

The charging station accommodates up to four central units, though the central units (motors) are not included with the charging station. When a central unit is inserted, it is checked by an intelligent control system. A display then indicates the battery's charge status.










#### Cordless single use shaver system for arthroscopy

### Disposable shaver

- No preparation or sterilization required
- Ready for immediate use
- Integrated tube connection
- Minimum order quantity: 5 pieces



- Aggressive meniscus resection, fast removal of soft tissue during shoulder surgeries.
- Suitable for a combination of soft tissue and bone resection.

Ø	Item No.
3,5mm	SFR-103L
5,0mm	SFR-104L



Versatile aggressive cutting blade: Ideal for extensive tissue removal.

Ø	Item No.
2,9mm	SFR-102N
5.0mm	SFR-104N



- Oval Burr: Fast bone resection.
- Suitable for subacromial decompression.

Ø	Item No.
4,5mm	SFR-104E



 Serrated inner blade and smooth outer blade for effective resection of meniscus, cartilage, and synovium.

Ø	Item No.	
3,5mm	SFR-103P	



Ball Burr: Ideal for notchplasty and osteochondral cleaning.

Ø	Item No.
4,5mm	SFR-104C





#### **Camera Covers**

DE-10-25 **Camera Covers** 

> disposable, 1 pcs, polyethylen, 1 3x250 cm, minimal amount 50 pcs





#### **Cleaning accessories**

**DE-5-17** Cleaning brush, small, 17 cm

for tubes, ID:2-4 mm, autoclavable

DE-10-300 **Disinfection solution IP** 

> for flexible and rigid endoscopes. manual and automated use, 1-3% dilution, 5 litres



Soak basin DE-10-425

> for chemical solution disinfection, cons. of basin (perforated) tray & cover for basin, size: 42x16x9,5 cm







#### Video cart

Completely mounted

- √ 50cm unit size
- ✓ Camera fixation knob
- ✓ Drawer
- ✓ Ask for additional accessories

WV-54043 Video cart, H: 143 cm

size: 62x63 cm, completely mounted, basic version,

color: grey, partly blue, \*\*drawer inbuilt

\*\*fixation knob for video camera head

\*\*specially rust protected max. weight: 150kg/330lbs

(Gastroscope bracket, fixation system. devices and monitor are separate accessory)

WV-53022 Video cart

as above but height 110 cm (without pictures)



WV-54043

WV-45111

Video cart, H: 112 cm,

size:54x58cm,

completely mounted, basic version,

drawer inbuilt with handle max.weight:50kg/110lbs max. unit size: 45cm color:grey, partly blue,



#### **Accessories for Video cart**

WV-99619 **Fixation system for infusion bottles** 

mounting set

WV-99618G **Gastroscope mount for Dr. Fritz TG series** 

inductive charging, with telescope rod

WV-01601 Fixation system for CO<sup>2</sup> bottles

> mounting set, comp.to video carts for 5-10 litre bottles, max OD: 14,5 c

WV-02851 Fixation system for CO<sup>2</sup> bottles

> mounting set, comp to video carts for big bottles, max OD: 20,5 cm








## Dogs, Cats, Birds & Exotic Species

**Endoscopic equipment** 



# Endoscopy

in Pets, Birds and Exotic Species





## Dogs, Cats, Birds & Exotic Species

**Endoscopic equipment** 

### **Video Compact System 7**

More information on page 10

- √ Wireless video transmission
- √ Independent of video tower and light source
- √ Without cable
- √ Easy handling & quick cleaning
- √ Integrated LED light

Chip-in-the-Tip-Technologie robust and unbreakable



General use for small Animal consulting rooms

- Rhinoscopy
- Laryngoscopy
- Dental endoscopy
- Otoscopy
- Vaginoscopy
- Dermoscopy
- Best compliance with the owner

#### Set 98-VCS7-SE

- VCS7, Ø 3,4mm, WL: 6cm, working channel Ø1,2mm
- 10" monitor receiver
- Application sheaths for ears and oral examination in pets
- Optical grasping forceps
- Flexible grasping forceps Ø1mm, WL: 16cm
- · Ear curette
- Charging station, wall mount
- Antifog, cleaning brushes
- Stainless steel stand for silicone sleeves including charging station
- Inflation pump for dental and gyn







## Dogs, Cats, Birds & Exotic Species

**Endoscopic equipment** 

#### **Compact Rhinoscope/Vaginoscope**

More information on page

#### Set (TA20-3SET)

- Compact Rhinoscope/Vaginoscope
- 2 spare batteries
- Charging station
- Charger
- Biopsy forceps
- 3 way stopcock
- Catheter, 1,0mm
- Cleaning brush
- Protection sheath NEW: enables circulation flushing
- 10" monitor receiver

Cordless generation with exchangeable battery, in-built LED's, video camera, HD sender for video transmission & 1,5" control display.



#### **Video Bronchoscope**

More information on page

#### Set TXBH4-54V

- Bronchoscope TXBH4-54
- Charger
- Cleaning brush
- Density tester
- Grasping forceps
- Storage case
- 7" Monitor receiver





Item No.	Length	Ø	Working channel	Transmission
TXBH2-36V	36cm	2,9mm	No working channel	wireless
TXBH4-54V	54cm	3,9mm	1,2mm	wireless

