

Equine



Video- Gastroscope

Professional gastroscopes - connection for insufflation, irrigation and suction



WIRELESS Video Gastroscope

- NextGen4 -





Video- Gastroscope

Professional gastroscope - connection for insufflation, irrigation and suction



- ✓ Integrated sender for wireless video transmission
- ✓ Integrated video camera with LED light
- ✓ Integrated DVR
- ✓ Ergonomic body
- ✓ Easy handling, quick cleaning and disinfection


New Generation - NextGen4

- Integrated control display
- Integrated SD card slot
- Fully soakable
- Truly mobile

Item No.	Length	Ø	Working channel	Transmission
TG326-12	320cm	12mm	3,2mm	Wireless

Set TG326-12G

- Gastroscope TG326-12
 - Charger
 - Cleaning adapters
 - Cleaning brush
 - Density tester
 - Biopsy forceps
 - 10" Receiver-monitor, recorder
 - Air/water pump
 - Transportation case
- Optional: wall mount



Wireless HD transmission





■ Set TG326-12G

Gastroscope
OD: 12mm
WL: 3,20m
(TG326-12)



Charger
magnetic
(TV27-NET)




Cleaning adapters
(TW132-24S)



Cleaning brush
(F28-3008)



Biopsy forceps
(B23-3601)



Density tester
(EF-23)



Air/water pump
(TW4-250SE)



10" Receiver-monitor,
recorder
(TXM10-58R)



Case for NextGen4
Size: 56x41x11cm
(K-5641-10)





Additional equipment


Video Gastroscope

Air/water pump

For lens cleaning and inflation

TW4-250SE

- Air/water pump (H150-1-1)
- Air/water twin tube (TW4-133)
- Fixation for bottle (TW4-261H)
- Bottle 250ml (TW4-250)




Suction-/irrigation pump

For lens cleaning, inflation and suction.
For small and large animal practice.

H140-3S

- Suction & irrigation pump (H140-3)
- 1 li. suction bottle (H160-AB-43)
- 1 li. flushing bottle (H160-AB-7)
- Tubes & transformer
- Sterile filter (AR-10-922)




Wall mount

Wall mount for inductive charging
- always ready for use, perfectly
stored and optimally drying.

TGX-25S

- Wall mount
- Power supply
- Mounting material





Cleaning accessories

Video Bronchoscope

■ For rigid & flexible endoscopes / instruments

DE-10-702 Dr. Fritz enzymatic cleaning FLEX 1,5 %,
1 x 1 l, 7 enzymes, 0,5%
specially for flexible endoscopes
for daily manual cleaning



DE-10-300 „IP“ Disinfection solution
for flexible and rigid endoscopes, manual and
automated use, 1-3% dilution, 5 litres



■ Option for cleaning of flexible endoscopes in the field



- ✓ Close the lid and ready for transport
- ✓ Fast and easy pre-cleaning after the last examination
- ✓ For semi- or fully automatic cleaning
- ✓ For all flexible endoscope types
- ✓ Usefull in drawer cabinets
- ✓ With airing



Stackable



DE-17-44S Transportation & storage basin
for flexible endoscopes, 50x44x16 cm,
set with cover

consisting of:

- Basin (DE-17-44)
- Cover transparent (DE-17-41)