


Equine



# Compact Video Vaginoscope

Professional video vaginoscope, wireless, for diagnostic use and treatment



# **WIRELESS** Compact Video Vaginoscope

worldwide, unique and up to now  
unmatched easy in practical use



Control display



Battery replaceable

97 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com

**Dr. Fritz**<sup>®</sup>  
Wireless ■ Endoscopes



# Compact Video Vaginoscope

Professional video vaginoscope, wireless, for diagnostic use and treatment

## For practitioners and universities



- ✓ Battery replaceable
- ✓ Control display on the endoscope
- ✓ Integrated sender 5,8 GHz
- ✓ Integrated LED light
- ✓ Easy handling & quick cleaning
- ✓ Fully immersible

### Set TQ52-28PF

Compact Video Vaginoscope & speculum for equine. Cordless generation with exchangeable battery, built-in camera with 5,8 GHz radio transmission, LED light, 1,5" display. Both completely soakable when vaginoscope and battery are separated.

- Compact Video Vaginoscope
- Obturator, blunt
- 2 spare batteries
- Charging station
- Charger
- Reducer adapter
- 21" monitor receiver

## Accessories to be used as transvaginal laparoscope

### TQ52-65

**Fornix perforator / cross trocar**, 26x570mm. Compatible with TQ52-28PF.

- Trocar & blade guide (TQ52-65-1/2)
- Spare blades, 12 pcs./pack (TQ52-66KL)
- Spare part set, hexagonal wrench, spring, protection caps (TQ52-65ER)



## Spare parts

### Reducer adapter

for TQ52-28, Luer (fem), cap w.5-10mm (234-211), opening and cover cap (234-600) (TQ52-E)

### Scalpel blades

No.20 for trocar TQ52-65, individually packed, 12pcs/pack (TQ52-66KL)



## Individual parts TQ52 series

Compact Video Vaginoscope  
OD: 28mm , WL: 52cm  
fully immersible  
(TQ52-28)



Charger  
(97-FU14EU)

Charging station  
(TVG-4)

Spare Li-Battery  
(TVG-36)



Obturator, blunt  
(TQ52-281)



Reducer adapter  
luer (fem), cap w.5-10mm,  
opening and cover cap  
(TQ52-E)



Scalpel blades, 12pcs.  
(TQ52-66KL)



Fornix perforator/ cross trocar  
(SET: TQ52-65)



Case,  
Case lid can be used as a  
disinfection tray  
Dimensions: 84x27x11cm  
(Y-K-8427-11-1 /  
Y-TQ52-SC)



21" Monitor receiver  
incl. recorder & accu  
(M-LCD-21SK)

