


Small Animal



Compact Rhinoscope-/Vaginoscope

Classic endoscopic equipment for pets- and small animals



WIRELESS Compact Rhinoscope/Vaginoscope

One of the smallest units with
working channel worldwide

15 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com

Dr. Fritz[®]
Wireless ■ Endoscopes



Compact Rhinoscope/Vaginoscope

Classic endoscopic equipment for pets and small animals



- ✓ Battery replaceable
- ✓ Control display on the endoscope
- ✓ Integrated sender 5,8 GHz
- ✓ Integrated LED light
- ✓ Easy handling & quick cleaning
- ✓ Fully immersible



Control display



Battery replaceable

■ Set (TA20-3SET)

Cordless generation with exchangeable battery, in-built LED's, video camera, HD sender for video transmission & 1,5" control display.

- Compact Rhinoscope/Vaginoscope
- 2 spare batteries
- Charging station
- Charger
- Biopsy forceps
- 3 way stopcock
- Catheter, 1,0mm
- Cleaning brush
- Protection sheath
- NEW: enables circulation flushing
- 10" monitor receiver

▷ Several applications

see page 13

Wireless video transmission





Compact Rhinoscope/Vaginoscope

Classic endoscopic equipment for pets and small animals

■ Set (TA20-3SET)


Compact Rhinoscope-/Vaginoscope
OD: 3,4x2,6mm, 8,5 Charr
WL: 20cm
Working channel: 1,2mm
Luer connection, fully immersible



Charger
(97-FU14EU)

Charging station
(TVG-4)

Spare Li-Battery
(TVG-36)



Biopsy forceps
OD: 1,0 mm
WL: 35cm
(B10-351)



Flushing & protection cannula
4mm
enables 2 way suction & irrigation
(TV20-3SP)



Cleaning brush



3 way stopcock
(Y-LL-1273H)



Catheter
OD: 1,0 mm
WL: 100cm
(F13-102A)



10" Monitor receiver
incl. recorder & accu
(TXM10-58R)





Additional equipment


Compact Video Rhinoscope/Vaginoscope

Suction-/irrigation pump

H140-3S

- Suction & irrigation pump (H140-3)
- 1 li. suction bottle (H160-AB-43)
- 1 li. flushing bottle (H160-AB-7)
- Tubes & transformer
- Sterile filter (AR-10-922)

For lens cleaning, inflation and suction.
For small and large animal practice.



NEW! Post rhino video mirror

TA20-3P18

OD: 18mm, WL: 20cm
With soft palate hook

Until now examination of retro pharynx was only possible for clinics and hospitals with appropriate endoscopic equipment. The newly developed video mirror system allows the general practitioners to perform retro pharyngoscopy. Grasping of foreign bodies or biopsies under visual control are possible.






Post rhinoscopy LED mirror system

for cats and dogs



Until now examination of retro pharynx was only possible for clinics and hospitals with appropriate endoscopic equipment. The newly developed deflecting mirror system allows the average practice to perform retro pharyngoscopy. Grasping of foreign bodies or taking biopsies under visual control are possible.

Choanae mirror with LED light, independent of endoscopes



E-24-LE-SET Post rhinoscopy LED set consisting of:

- E-24-LE Light carrier
- E-24-21 Choanae mirror
- E-22-2 Choanae hook
- E-22-71 Antifog

Zubehör



E-22-2 Epiglot double hook
16 cm, for cats and middle sized dogs



E-22-71 Antifog fluid
for endoscopes and mirrors, 50 ml

K-3325-61 Special transportation case
33x25x6 cm (without picture)