





Dogs, Cats, Birds & Exotic Species

Endoscopic equipment



Endoscopy

in Pets, Birds and Exotic Species





Dogs, Cats, Birds & Exotic Species

Endoscopic equipment

Video Compact System 7

More information on page 10

- ✓ Wireless video transmission
- ✓ Independent of video tower and light source
- ✓ Without cable
- ✓ Easy handling & quick cleaning
- ✓ Integrated LED light

NEW ⇒ Integrated working channel Ø 1,2mm
Ø 3,4mm & WL: 6cm

NEW ⇒ Chip-in-the-Tip-Technologie
robust and unbreakable



General use for small Animal consulting rooms

- Rhinoscopy
- Laryngoscopy
- Dental endoscopy
- Otoscopy
- Vaginoscopy
- Dermoscopy
- Best compliance with the owner

Set 98-VCS7-SE

- VCS7, Ø 3,4mm, WL: 6cm, working channel Ø1,2mm
- 10" monitor receiver
- Application sheaths for ears and oral examination in pets
- Optical grasping forceps
- Flexible grasping forceps Ø1mm, WL: 16cm
- Ear curette
- Charging station, wall mount
- Antifog, cleaning brushes
- Stainless steel stand for silicone sleeves including charging station
- Inflation pump for dental and gyn





Dogs, Cats, Birds & Exotic Species

Endoscopic equipment

Compact Rhinoscope/Vaginoscope

More information on page

Set (TA20-3SET)

- Compact Rhinoscope/Vaginoscope
- 2 spare batteries
- Charging station
- Charger
- Biopsy forceps
- 3 way stopcock
- Catheter, 1,0mm
- Cleaning brush
- Protection sheath
NEW: enables circulation flushing
- 10" monitor receiver

Cordless generation with exchangeable battery, in-built LED's, video camera, HD sender for video transmission & 1,5" control display.



Video Bronchoscope

More information on page

Set TXBH4-54V

- Bronchoscope TXBH4-54
- Charger
- Cleaning brush
- Density tester
- Grasping forceps
- Storage case
- 7" Monitor receiver



Item No.	Length	ø	Working channel	Transmission
TXBH2-36V	36cm	2,9mm	No working channel	wireless
TXBH4-54V	54cm	3,9mm	1,2mm	wireless