

Equine



Compact Video Vaginoscope

Professional video vaginoscope, wireless, for diagnostic use and treatment



WIRELESS Compact Video Vaginoscope

worldwide, unique and up to now
unmatched easy in practical use



Control display



Battery replaceable

99 2025

+49 (0)7777 / 882300 ■ dr.fritz@dr-fritz.de ■ www.dr-fritz.com

Dr. Fritz[®]
Wireless Endoscopes



Compact Video Vaginoscope

Professional video vaginoscope, wireless, for diagnostic use and treatment

For practitioners and universities



- ✓ Battery replaceable
- ✓ Control display on the endoscope
- ✓ Integrated sender 5,8 GHz
- ✓ Integrated LED light
- ✓ Easy handling & quick cleaning
- ✓ Fully immersible

Set TQ52-28PF

Compact Video Vaginoscope & speculum for equine. Cordless generation with exchangeable battery, built-in camera with 5,8 GHz radio transmission, LED light, 1,5" display. Both completely soakable when vaginoscope and battery are separated.

- Compact Video Vaginoscope
- Obturator, blunt
- 2 spare batteries
- Charging station
- Charger
- Reducer adapter
- 21" monitor receiver

Accessories to be used as transvaginal laparoscope

TQ52-65

Fornix perforator / cross trocar, 26x570mm. Compatible with TQ52-28PF.

- Trocar & blade guide (TQ52-65-1/2)
- Spare blades, 12 pcs./pack (TQ52-66KL)
- Spare part set, hexagonal wrench, spring, protection caps (TQ52-65ER)



Spare parts

Reducer adapter

for TQ52-28, Luer (fem), cap w.5-10mm (234-211), opening and cover cap (234-600) (TQ52-E)

Scalpel blades

No.20 for trocar TQ52-65, individually packed, 12pcs/pack (TQ52-66KL)



Individual parts TQ52 series

Compact Video Vaginoscope
OD: 28mm , WL: 52cm
fully immersible
(TQ52-28)



Charger
(97-FU14EU)

Charging station
(TVG-4)

Spare Li-Battery
(TVG-36)



Obturator, blunt
(TQ52-281)



Reducer adapter
luer (fem), cap w.5-10mm,
opening and cover cap
(TQ52-E)



Scalpel blades, 12pcs.
(TQ52-66KL)



Fornix perforator/ cross trocar
(SET: TQ52-65)



Case,
Case lid can be used as a
disinfection tray
Dimensions: 84x27x11cm
(Y-K-8427-11-1 /
Y-TQ52-SC)



21" Monitor receiver
incl. recorder & accu
(M-LCD-21SK)

